



**NEW.
1 MINUTE
KILL TIME FOR
SARS-CoV-2**

NAVIGATOR

Dilution Control System


Why It's A Better System

- Concentrate and RTU labels are numbered and color coded for safety and to facilitate training
- Closed system - Spill proof container
- New Metered Tip with each concentrate provides proper dilution each time
- Each concentrate has an easy carry handle
- Recyclable container

Why It's A Better Cleaner

- Effective against SARS-Related Coronavirus 2² (SARS-CoV-2) in 1 minute
- 3x Renegade cleans, disinfects and deodorizes in one step
- Disinfectant designed for Hospital Use
- Effective as a Virucidal and against a broad spectrum of bacteria including Human Immunodeficiency Virus Type 1 (HIV-1) associated with AIDS or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV)
- Provides Norovirus efficacy in fast 5 minute contact time
- Fungicidal, inhibits the growth of mold and mildew, and their odors, when used as directed
- Helps reduce cross-contamination and healthcare acquired infections (HAIs)
- Neutral pH won't dull finished floors

How It Cleans



FOR USE IN A SPRAY BOTTLE: Prepare a fresh solution daily. Dispense diluted product into properly labeled RTU dry container. For heavily soiled areas a preliminary cleaning is required.



Spray 6-8 inches from a hard, non-porous surfaces, thoroughly wetting the surface with a sprayer. Treated surfaces must remain wet for contact time indicated on the label.



Allow to air dry. Do not breathe spray. Rinse food contact surfaces such as countertops, appliances, tables with potable water before reuse.



FOR USE IN AN AUTOMATIC SCRUBBER OR MOP BUCKET: Dispense diluted product into mop bucket or automatic scrubber. Apply product on surfaces to be cleaned. Treated surfaces must remain wet for contact time indicated on the label.

FOR USE IN MICROFIBER CHARGING SYSTEMS
Dispense diluted product into microfiber mop bucket according to chart below. Charge microfiber until completely saturated (15-20 minutes). Treated surfaces with charged microfiber cloth should remain wet for contact time indicated on label.

Type	Qty	Diluted Solution
16"x16" microfiber cloths	15	32 oz
	24	64 oz
	40	96 oz
	50	1 gal
9" wall wash pads	20	96 oz
	25	1 gal
	30	1.25 gal
	35	1.3 gal
18" microfiber pads	10	96 oz
	15	1.3 gal
	20	1.5 gal
	25	1.75 gal

3x RENEGADE DAILY

ONE-STEP DISINFECTANT
KILLS GERMS IN 3 MINUTES - LABOR SAVING

EPA Registered Disinfectant

EPA Reg. No. 6836-349-12120
EPA EST. No. 8325-PA-01

Where It Cleans

- Floors and Walls
- Counter Tops
- Tables and Chairs
- Bathroom Fixtures
- Metal Surfaces
- Showers
- Noncritical Medical Devices
- Bed Rails
- Other Hard Nonporous Surfaces

Technical Data

Packaging (concentrate) 2L, 3L PDC, 1G
Product Number 13215,48333,48407
Dilution: 1:128
Color: Yellow
Odor: Lemon
pH: 12-13.5
RTU pH: 9.5
Foam: Moderate
Solubility: Complete with water
Detergency: Excellent
Rinsability: Excellent
Biodegradability: Complete
Storage/Stability: 1 year
Freeze/Thaw Stability: Keep from freezing
RTU Quart Yield: 270 Quarts per bottle
Gallon Yield: 67 Gallons per bottle
PDC RTU Quart Yield 405 Quarts per bottle
PDC Gallon Yield 101 Gallons per bottle

Navigator Dispensers

- #23200 Navigator DCS 4-Product Dispenser
- #23201 Navigator DCS 1 Product Dispenser, bottle-fill
- #23202 Navigator DCS 1 Product Dispenser bucket-fill
- #23203 Navigator DCS Portable Dispenser
- #23204 Navigator DCS Portable w/Basic Spray/Foaming Dispenser
- #23205 Navigator DCS Foaming Dispenser w/Ext. Hose
- #23209 6 ft. Water Supply Hose, Black
- #23215 Navigator MPD Wall Mount Station - Bucket & Bottle fill

